



# TECHNICAL DATASHEET

IM-30-04/1

Edition: 13/05/13

Page: 1 of 2

## PRODUCT NAME:

NEOPRENE CROSSOVER ELASTIC ANKLE SUPPORT- THERMO-MED®

## REFERENCE:

4401 (Blue) – 5401 (Beige)

## SIZES:

SIZES	PERIMETER CM	
XS/1	22-24	
S/2	24-26	
M/3	26-28	
L/4	28-31	
XL/5	31-34	

## MANUFACTURER:

ORLIMAN, S.L.U. CIF: B-96122510  
C/ AUSIAS MARCH Nº 3. POLIGONO INDUSTRIAL LA POBLA – L'ELIANA.  
46185 LA POBLA DE VALLBONA – VALENCIA – SPAIN  
SPANISH MEDICAL DEVICE OPERATION LICENCE NUMBER: 20-PS

## PRODUCT CLASSIFICATION (Council Directive 93/42/EEC (as amended by Council Directive 2007/47/EC))

CLASS: I – (Rule 1)

Manufacturing type: Mass-Produced.

## PRODUCT CHARACTERISTICS:

3 mm neoprene ankle support with trimming parts in the upper part for better adjustment

### Colours:

- Blue: 4401
- Beige: 5401

## INDICATIONS:

- Slight sprains.
- Tendonitis.
- Post-operatives.
- Relapses.

## MATERIALS:

- Components: Contain LATEX, natural rubber.
- Fabric: Inside: Polyamide 100%; Intermediate: Polychloroprene 100%; Outside: Polyamide 100%.
- Strap: Elastic neoprene.

## PRODUCT PICTURE:



Difarza, s.l. CIF B96386875 Tel.: 963363215 Fax: 963671551  
email: [disalud@ono.com](mailto:disalud@ono.com)  
[www.disalud.com](http://www.disalud.com)

C/Río Ulla, 5 (entrada por Juan Verdeguer, 25) 46024 Valencia.  
Empresa colaboradora de servicios médicos y de prevención de riesgos.  
Distribución de productos sanitarios, clínicos y ortopédicos desde 1980.  
Prevención, Seguridad y Salud Laboral.



1112&10#%&2.\*3



1112#(3#.!#&10#%&2.\*3



!''#%&' (%) \*+, -#\$+, .!#



/ &!0#%&



## TECHNICAL DATASHEET

IM-30-04/1

Edition: 13/05/13

Page: 2 of 2

### RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.



### ADDITIONAL INFORMATION:

Guarantee: 6 month.

CE

Difarza, s.l. CIF B96386875 Tel.: 963363215 Fax: 963671551  
email: [disalud@ono.com](mailto:disalud@ono.com)  
[www.disalud.com](http://www.disalud.com)

C/Río Ulla, 5 (entrada por Juan Verdeguer, 25) 46024 Valencia.  
Empresa colaboradora de servicios médicos y de prevención de riesgos.  
Distribución de productos sanitarios, clínicos y ortopédicos desde 1980.  
Prevención, Seguridad y Salud Laboral.



[www.disalud.com](http://www.disalud.com)



[www.farmaciadisalud.com](http://www.farmaciadisalud.com)



DiSaludGrupoenValencia



@disalud